

# WHEELCHAIRS

## BASIC MANUAL WHEELCHAIRS



## NEW CLASSIC

### Manual, foldable, steel wheelchair

Standard wheelchair on a cross frame made of chrome-molybdenum steel. It has lifting armrests and footrests that are adjustable in length, rotating and removable. The upholstery is made of non-flammable material in black. The frame is powder coated in blue. The wheelchair is equipped with black, specially profiled push rims, that make it easier for user to move around independently. The New Classic can be attached easily to a caterpillar stairclimber.

Ask for options:

1. Rear wheels with quick release axles.
2. Hip safety belt
3. Hand brakes for an attendant



Front wheels size	Rear wheels size
200 x 33 mm (PU) – puncture proof	24" x 1 3/8" pneumatic Can be replaced with PU, puncture proof tires

Seat width	Seat depth	Backrest height	Overall width	Seat to floor height	Seat to footplate height	Overall weight	Max weight of user
40, 45, 50 cm	40 cm	44 cm	59, 64, 69 cm	48 cm	35-46 cm	17,35 kg	120 kg