

WHEELCHAIRS

BASIC MANUAL WHEELCHAIRS



FLIPPER

Manual, foldable, aluminum, light weight wheelchair with quick release axles

Aluminum wheelchair on a cross, foldable frame, equipped with quick release mechanism for rear wheels. It has height adjustable, lifting and detachable armrests. Its footrests are adjustable in length, rotating and removable. Flipper is equipped with an assembly plates which allow to adjust position of rear wheels. It is recommended for patients after leg amputations. Flipper represents known European brand.

Flipper can be attached easily to a caterpillar stairclimber.



Front wheels size	Rear wheels size
200 x 37 mm (PU) – puncture proof	24 x 1 3/8 pneumatic Can be replaced with PU, puncture proof tires

Seat width	Seat depth	Backrest height	Overall width	Seat to floor height	Seat to footplate height	Overall weight)	Max weight of user
40, 44, 48, 51 cm	44 cm	45 cm	61,5; 65,5; 69,5; 72,5; cm	49 cm	34-48 cm	16,2 kg	136 kg