

# WHEELCHAIRS

## BASIC MANUAL WHEELCHAIRS



## DOLPHIN

**Manual, foldable, aluminum, light weight wheelchair with quick release axles**

Aluminum wheelchair on a cross, foldable frame, equipped with quick release mechanism for rear wheels. It has height adjustable, lifting and detachable armrests. Its footrests are adjustable in length, rotating and removable. As a standard feature the wheelchair is equipped with height adjustable push handles for an attendant and with seat cushion. It is possible to adjust the seat depth.. Dolphin is equipped with an assembly plates which allow to adjust position of rear wheels and its front forks include angle adjustment.. It is recommended for patients after leg amputations. Dolphin represents known European brand.



Front wheels size	Rear wheels size
200 x 45 mm (PU) – puncture proof	24 x 1 3/8 (PU) – puncture proof

Ask for optional equipment: drum brakes with hand levers for an attendant  
Ask for optional version with max load increased up to 200 kg in 54 cm & 57 cm seat widths.

Seat width	Seat depth	Backrest height	Overall width	Seat to floor height	Seat to footplate height	Overall weight)	Max weight of user
40, 44, 48, 51 cm	42-48 cm	44 cm	60,3; 64,3; 68,3; 71,3 cm	45-52,5 cm	34-48 cm	16,75 kg	150 kg