

WHEELCHAIRS

MANUAL WHEELCHAIRS



Optional equipment

Anti tip



Lift up footrest



Drum brakes



Table



BARRACUDA

Special wheelchair, with back and head stabilization and quick release axels on rear wheels

Aluminum wheelchair with reclining backrest and head support on a cross frame for full day care. The wheelchair is equipped with backrest with angle regulation (30°) and possibility of individual adjustment of seat depth, seat width, height of front and back of the wheelchair. Height of armrests, legrests and attendant handles can be regulated. Thank to rear wheels adapters the center of gravity of Barracuda can be easily changed and adapted for patients after limbs amputation.

Standard features:

- Seat depth adjustment (42 – 48 cm).
- 6 position assembly plates of rear wheels and 4 position adapter of front wheel assembly (rear and front wheels equipped with quick release system)
- Fluent backrest angle regulation (30°)
- Height and angle regulated headrest (360°)
- Anti tip
- Height regulated and swing away legrests
- Swing away, detachable and height regulated armrests
- Push handles with independent height adjustment
- Seat cushion
- Adjustable backrest upholstery tension

Legrest length	Armrest height	Backrest recline
36,5 – 52 cm	17,5-26 cm	90° – 120°

Front wheels size	Rear wheels size
200 x 45 mm (PU) – puncture proof	24" x 1 3/8" (PU) – puncture proof

Seat width	Seat depth	Backrest height	total width	Seat to floor height	Total length with legrest	Weight of wheelchair depending on size	Max load
40, 44, 48, 50 cm	42-48 cm	38 cm	60,3; 64,3; 68,3; 71,3 cm	45-52,5 cm	99,5-105,5 cm	20-25 kg	120 kg