

WHEELCHAIRS

POWER SCOOTERS FOR THE DISABLED



Go Go Elite Traveller

Power scooter for the disabled

GO GO Elite Traveller scooter is the easiest solution for your personal independence and comfort of your everyday travelling. GO GO scooters were specially designed for elderly people and for persons who are partly disabled, so they can easily fulfill their daily routines.

GO GO is a compact and agile scooter, very easy to maneuver, even in narrow spaces. Overall weight of the device, without batteries is below 40 kg. The scooter, almost without folding fits in to trunks of many cars which makes it very handy in transportation or travelling (e.g. for holidays, to a shopping center, to your recreational cottage, to a park or anywhere where your GO GO may be needed). Spinal or joints diseases, muscle pains, partial paresis, cardiac or venous problems won't be a holdback anymore when going out for a walk or shopping. GO GO delivers independence of daily moving around. Remember, GO GO can be taken anywhere and using it is pure fun and freedom. GO GO can be delivered in one of two versions: as a 4-wheeler or 3-wheeler. The most advanced model (GO GO Elite Traveller PLUS) has higher max. load capacity: 147 kg and loop shaped steering handle for easier maneuvering.

TECHNICAL DATA	Go Go Traveller	Go Go Traveller Plus
Batteries	2 x 12 V	2 x 12 V, 17 Ah
Max. speed	7 km/h	7,5 km/h
Range	16 km	17,5 km
Max. inclination of slopes	8°	8°
Width	100,3 cm	102 cm
Weight without batteries	34 kg	39,5 kg
Weight of batteries	9,5 kg	13 kg
Max. weight of a user	125 kg	147 kg
Turning radius	112 cm	116 cm
Front wheels size	8" – PU, puncture proof	8" – PU, puncture proof
Rear wheels size	9" – PU, puncture proof	9" – PU, puncture proof
Ground clearance	7,62 cm	7,62 cm