

WALKERS & ROLLATORS

ROLLATORS



JAGUAR ROLLATOR

Jaguar, is modernly designed rollator which is one of its kind..

It is a very special product dedicated for users who like its unique, colorful design. At the same time the rollator is fully loaded with equipment, maneuverable and compact.

Jaguar has aluminum frame, bag and tray. Max. user's weight is 100 kg and max load of the bag is 5kg. The rollator is very easy to fold for storage or transportation and it weighs only 6,6 kg. Main purpose of Jaguar is assistance in standing up, walking and seating down for disabled or elderly persons. It is Push handles are height adjustable for optimal user's fit.

TECHNICAL DATA

Frame material	aluminum
Push handles height (range)	83,5-94,5 cm
Width	54 cm
Length (folded)	25 cm
Length (opened)	66 cm
Wheels diameter	7"
Max. weight of user	100 kg
Max. load of the bag	5 kg
Weight	6,6 kg
Brakes	2 maneuvering brakes with parking lock function