

# WALKERS & ROLLATORS

## ROLLATORS



## LEOPARD ROLLATOR

LEOPARD is modern, beautifully designed 4 wheel rollator for the most demanding customers. It is made of aluminum and light (weights approx. 7 kg). The device assists disabled persons during standing, seating or walking. It is handy, agile and very functional. Max. weight of a user is 150 kg and bag can be loaded with belongings weighing up to 5 kg. Thanks to its cross bar frame design Leopard is easy to fold for storage or transportation.

Leopard is equipped with soft seat which allows user to rest during longer walks in the park or during shopping. It has height adjustable push handles for optimal user's fit.

### TECHNICAL DATA

Frame material	aluminum
Push handles height (range)	79-99 cm
Width	60,8 cm
Width (folded)	22,5 cm
Length	66 cm
Seat to floor distance	61 cm
Seat width	46 cm
Wheels size	202 x 35 mm
Max. weight of user	150 kg
Max. load of the bag	5 kg
Weight	7,3 kg
Brakes	2 maneuvering brakes with parking lock function