

WALKERS & ROLLATORS

HIGH SUPPORT WALKERS



TIGER ROLLATOR

Modern and beautifully designed rollator with arm supports for demanding patients. Light aluminum construction allows to assist patients with disabilities while standing or walking. Rollator is very maneuverable, handy and functional. Its max load is 150 kg and the bag can carry load up to 5 kg. Cross folding allows easy folding for transportation and storage.

Tiger is equipped with softly padded height adjustable arm supports with easy reachable brakes. Soft and comfortable seat allows patients to rest while longer walks.

TECHNICAL DATA

Frame material	aluminium
Push handles height (range)	95 – 120 cm
Width	62 cm
Width (folded)	22,5 cm
Length	66 cm
Seat to floor distance	62 cm
Seat width	46 cm
Wheels size	200 x 35 mm
Max. weight of user	150 kg
Max. load of the bag	5 kg
Weight	9,5 kg
Brakes	2 maneuvering brakes with parking lock function