

# CANES & CRUTCHES

## CANES



### CANES WITH HEIGHT ADJUSTMENT FASHION

Aluminum, fashion canes facilitate walking for elderly or disabled persons. They are equipped with height adjustment function and with handgrip made of maple wood. Canes are available in 4 styles: scotch, flowers, chessboard, milk chocolate.



| Height adjustment (range) | Height adjustment steps | Weight | Internal diameter of rubber foot | Diameter of rubber foot's base | Max. weight of user |
|---------------------------|-------------------------|--------|----------------------------------|--------------------------------|---------------------|
| 71-94 cm                  | 2,5 cm                  | 0,4 kg | 19 mm                            | 38 mm                          | do 100 kg           |



### CANES WITH HEIGHT ADJUSTMENT

Aluminum canes facilitate walking for elderly or disabled persons. They are equipped with height adjustment function and with handgrip made of wood. Canes are available in 3 colors: black, brown & graphite (marbled).



| Height adjustment (range) | Height adjustment steps | Weight | Internal diameter of rubber foot | Diameter of rubber foot's base | Max. weight of user |
|---------------------------|-------------------------|--------|----------------------------------|--------------------------------|---------------------|
| 71-94 cm                  | 2,5 cm                  | 0,4 kg | 19 mm                            | 38 mm                          | do 100 kg           |