

# CANES & CRUTCHES

## CANES



## FOLDABLE CANES

### 06.B.RPN.13

Aluminum, foldable cane facilitate walking for elderly or disabled persons. Thanks to its folding function it won't take much space in the bag or luggage, when there is no need to use it. The cane is equipped with height adjustment function and with handgrip made of wood. It is available in 2 colors: black & copper. The device is delivered with plastic case for storage.



\* ask for available colors

## SWAN CANE

Aluminum cane with height adjustment facilitate walking for elderly or disabled persons. It is equipped with specially designed soft hand grip which increases comfort of usage. Unique shape of the cane improves its balance. It is available in 2 colors:

- anodized gold
- anodized brown

Type	Height adjustment (range)	Height adjustment steps	Weight	Internal diameter of rubber foot	Diameter of rubber foot's base	Max. weight of user
foldable	82-92 cm	2,5 cm	0,35 kg	19 mm	38 mm	do 115 kg
swan	73-95 cm	2,5 cm	0,4 kg	19 mm	30 mm	do 100 kg



## WOODEN CANES

Wooden canes facilitate walking for elderly or disabled persons. They are available with one of 3 types of handgrip profiles: bowed, fritz & derby.



Height	Weight	Internal diameter of rubber foot	Diameter of rubber foot's base	Max. weight of user
93 cm	0,25 kg	19 mm	30 mm	do 100 kg