

# CANES & CRUTCHES

## RUBBER FEET FOR CANES & CRUTCHES



### RUBBER FEET FOR CRUTCHES

Rubber feet for crutches of different types. Available in following colors: black, grey & orange.



Internal diameter of rubber foot	External diameter of rubber foot's base	Rubber foot's height
19 mm (for elbow crutch foot) 22 mm (for underarm crutch foot & for hinged crutch foot)	38 mm (for elbow crutch foot) 42 mm (for underarm crutch foot & for hinged crutch foot)	50 mm (for elbow crutch foot) 55 mm (for underarm crutch foot & for hinged crutch foot)



### RUBBER FEET FOR CANES

Universal, rubber feet for canes of different types. Available in black.



Internal diameter of rubber foot	External diameter of rubber foot's base	Rubber foot's height
19 mm	38 mm	45 mm



### WINTER ANTI-SLIP TIPS FOR CRUTCHES & CANES

Universal anti-slip tip, which is easy to mount on all types of canes & crutches.

