

# BATHROOM EQUIPMENT

## SHOWER STOOLS & SEATS



### BATH BENCH WITH HANDLE

Universal bench for positioning on top of a bathtub, equipped with safety handle. It enables disabled persons to use bathtubs. It can be set for bathtubs with different widths.



Width	Depth	Height	Max. weight of user	Safety supports spacing (range)	Handle's height
70,5 cm	32 cm	18 cm	150 kg	41 - 64,5 cm	6 cm



### SHOWER STOOLS

Height adjustable shower stools. They are practical and nicely designed. They increase comfort of disabled or elderly persons during their daily hygienic activities. The triangular stool is equipped with soft seating pad for even better comfort. Thanks to its unique shape it takes very small space in bathroom or in shower cabin.



**SHOWER  
STOOL WITH  
BACKREST (SO)**



**SHOWER  
STOOL (SZ)**



**TRIANGULAR  
SHOWER  
STOOL**



**ROUND  
SHOWER  
STOOL**



**ROUND  
SHOWER  
STOOL WITH  
SWIVEL SEAT**

	Seat height adjustment range	Width	Depth	Max. weight of user
<b>STOOL (SO)</b>	40-58 cm	50 cm	30 cm	do 115 kg
<b>STOOL (SZ)</b>	40-58 cm	50 cm	30 cm	do 115 kg
<b>TRIANGULAR STOOL</b>	45,5-55,5 cm	40,4 cm	40,4 cm	150 kg
<b>ROUND STOOL</b>	40-50 cm	33 cm	33 cm	130 kg
<b>ROUND STOOL WITH SWIVEL SEAT</b>	34,5-52,5 cm	29,6 cm	29,6 cm	130 kg