

# BATHROOM EQUIPMENT

## RAISED TOILET SEATS



### RAISED TOILET SEAT WITH OR WITHOUT LID

This seat is easy to put on and to take off the toilet bowl at any moment, as it doesn't have to be attached permanently. It offers comfort and easiness of usage to all users. The seat is anatomically designed. With its 10 cm height it guarantees comfortable seating position and rubber, anti-slip pads stabilize its set up.



	Seat height	Width	Depth	Dimesions of the opening	Max. weight of user
Raised toilet seat with or without lid	10 cm	36 cm	40 cm	22 x 17,5 cm	do 150 kg



### RAISED TOILET SEAT WITH ARMRESTS (Fixed toilet seat raiser)

High toilet seat raiser with lid and with very stable toilet assembly. It is equipped with hinged, elevable armrests to enable an easy transfer on or off the toilet. Its height can be adjusted in following steps: 9/13/17 cm.



Seat height adjustment range	Width	Depth	Depth with armrests	Width between the top of armrests	Height of armrests above seat	Max. weight of user
9/13/17 cm	40 cm	49,5 cm	62 cm	53 cm	23 cm	do 120 kg