

STAIRCLIMBERS & RAMPS FOR THE DISABLED

RAMPS



DF/SF PORTABLE WHEELCHAIR RAMPS (foldable)

Our ramps are durable and light. Their unfolding or folding won't take more than 30 seconds. After folding ramps are not heavier or less handy than average suitcase. Thanks to an integrated handle ramps are easy to carry. Wide variety of available sizes allows to climb up different obstacles, from the lowest curbs up to 6 steps stairs.

FOLDING	LENGTH	MAX. LOAD	WEIGHT
single folding case SF	92 cm	300 kg	8 kg
single folding case SF	120 cm	300 kg	10,5 kg
single folding case SF	150 cm	300 kg	12,5 kg
double fold. case DF	180 cm	300 kg	16 kg
double fold. case DF	240 cm	300 kg	22 kg
double fold. case DF	300 cm	300 kg	27,2 kg