


STAIRCLIMBERS & RAMPS FOR THE DISABLED

RAMPS



THRESHOLD RUBBER RAMPS

Rubber ramps are easy solution to overcome thresholds in your living area. They are intentionally ribbed for anti-slip purpose and made of TPS rubber. Ramps are available with montage glue on the bottom or without. User can adjust them to individual needs by simply cutting them to desired length. Products are available in different sizes to eliminate obstacles with different heights.

Available ramp heights: 4, 6, 8, 10, 12, 16, 20, 24, 28, 32, 36, 40, 44, 48, 52, 56, 60 mm. Max. load is 300 kg.

Length of the riding surface depends on ramp height.

