

SPARE PARTS FOR WHEELCHAIRS

BATTERIES, INNERTUBES

BATTERIES



195 x 130 x 159 mm

197 x 165 x 170 mm

228 x 137 x 207 mm

259 x 168 x 208 mm

Gel battery
MWLG 33-12EV
12V 33Ah

Gel battery
MWLG 44-12EV
12V 44Ah

Gel battery
MWLG 55-12EV
12V 55Ah

Gel battery
MWLG 80-12EV
12V 80Ah

Maintenance free, MW Power gel batteries, of MWLG series are made with VRLA – GEL technology (gel electrolyte combined with PVC-SiO₂ separator including internal recombination of gases). They are designed for deep discharged working cycles such as: PV installations, electric vehicles, cleaning machines, power wheelchairs and golf carts, scooters for disabled persons, electric boats, camping power supply etc. This type of batteries can also be used for emergency power supply like: UPS, emergency illumination, automatic systems etc. Vitality of these batteries for 50% average discharge is up to 1000 working cycles. Service life is 12 – 15 years for work in temperatures 20 – 25° C.

INNERTUBES FOR WHEELCHAIRS



90/90 valve



90/30 valve



straight valve

Durable tubes for standard and active wheelchairs. The tubes are made of high quality rubber with high puncture resistance. Tubes are equipped with 1 of 3 types of car valves (as shown above) with black protection cap.

90/90 valve	90/30 valve	Straight valve
12 1/2x 2 1/4 (57x 203 mm)	6 x1 1/4 (150 x 30)	10 x 2 (52x 152 mm)
2.80/2.50-4	7 x 1 3/4 (175 x 45)	12 1/2 x 2 1/4 (57 x 203 mm)
3.00-4	8 x 1 1/4 (200 x 30)	22 x 1 (25 x 540 mm)
3.00-8	8 x 2 (200 x 50)	22 x 1 3/8 (37 x 498)
4.00-5		24 x 1 (25 x 540)
4.10/3.50-6		24 x 1 3/8 (37 x 540)
8x2 (200 x 50 mm)		3.00-8